



Model: HD-DIS001

Economic M Adult Diaper with Wetness Indicator

Suited waistline: 32 to 44-inch

Dimensions: 80 x 65cm

Pulp dimensions: 65 x 28cm

Weight: 95 ± 5g

SAP: 9 to 10g

Water absorbing capacity: 1,800ml

Packing: 10pcs/package



Model: HD-DIS002

Economic L Adult Diaper with Wetness Indicator

Suited waistline: 40 to 56-inch

Dimensions: 95 x 75cm

Pulp dimensions: 65 x 28cm

Weight: 105 ± 5g

SAP: 10 to 11g

Water absorbing capacity: 2,000ml

Packing: 10pcs/package



Model: HD-DIS003

Higher Absorbent M Adult Diapers with Frontal Tape

Suited waistline: 32 to 44-inch

Dimensions: 80 x 65cm

Pulp dimensions: 65 x 32cm

Weight: 120 ± 5g

SAP: 10 to 12g

Water absorbing capacity: 2,200ml

Front tape, stereoscopic leak guards

Packing: 10pcs/package



Model: HD-DIS004

Higher Absorbent L Adult Diapers with Frontal Tape

Suited waistline: 40 to 56-inch

Dimensions: 95 x 75cm

Pulp dimensions: 75 x 32cm

Weight: 140 ± 5g

SAP: 13 to 14g

Water absorbing capacity: 2,500ml

Front tape, stereoscopic leak guards

Packing: 10pcs/package



Model: HD-DIS005

Under Pad

Size 60 x 45cm, 43G/pc, absorbing capacity 500-600ml

Size 60 x 60cm, 60G/pc, absorbing capacity 700-800ml

Size 60 x 90cm, 100G/pc, absorbing capacity 1000-1100ml

Packing: 10pcs/package



Model: HD-DIS008 / HD-DIS009

Reinforced Endotracheal Tube

HD-DIS008, with cuff, I.D.5.0–9.5mm

HD-DIS009, without cuff, I.D. 3.0–4.5mm

Packing: sterile individual blister



Model: HD-DIS010 / HD-DIS011

Standard Endotracheal Tube

HD-DIS010, with cuff, I.D.3.0mm–10.0mm

HD-DIS011, without cuff, I.D.2.0mm–10.0mm,

Packing: sterile individual blister



HD-DIS012

Model: HD-DIS012

Laryngeal Mask Airway

Option: disposable or reusable

Size: 1.0mm 1.5mm 2.0mm 2.5mm 3.0mm 4.0mm 5.0mm



HD-DIS013

Model: HD-DIS013

Stomach tube

Size: 6 Fr-24 Fr

Option of tube: Transparent tube, Mist surface tube

X-Ray traceable available

Packing: Sterile individual polybag OR sterile blister

**Model: HD-DIS014**

Feeding tube

Size: 4 Fr–22 Fr

Option of tube: Transparent tube, Mist surface tube

X-Ray traceable available

Packing: Sterile individual polybag OR sterile blister

**Model: HD-DIS015**

Suction Tube

Size: 5 Fr–24 Fr

Option of tube: Transparent tube, Mist surface tube

Packing: Sterile individual polybag OR sterile blister


**Model: HD-DIS016**

Nelaton catheter

Size: 6 Fr–24 Fr

Option of tube: Transparent tube, Mist surface tube

Packing: Sterile individual polybag OR sterile blister

**Model: HD-DIS017**

Rectal Tube

Size: 24 Fr–36 Fr

Option of tube: Transparent tube, Mist surface tube


Packing: Sterile individual polybag OR sterile blister

**Model: HD-DIS018**

Redon drain catheter

Size: 6Fr-18Fr


Packing: sterile Individual polybag OR sterile blister

**Model: HD-DIS019**

Thoracic Duct Drainage

Size: 20Fr-36Fr

Packing: Sterile individual polybag OR sterile blister

**Model: HD-DIS021**

Penrose drainage tube

Size: 1/8"; 1/4"; 3/8"; 1/2"; 5/8"; 3/4"; 7/8"; 1"; 1.5"; 2"

Packing: sterile individual blister

**Model: HD-DIS022**

Pezzer Catheter

Size: 12Fr-34Fr, 400mm

Packing: sterile blister

**Model: HD-DIS023**

Latex T-drainage tube

Size: 8Fr-26Fr

Packing: sterile blister

**Model: HD-DIS024**

Latex foley catheter-3way standard, silicone coated

Size: 16Fr-26Fr, 30ml balloon

Available with the option of plastic valve or rubber valve

Packing: sterile individual blister

**Model: HD-DIS025**

Latex foley catheter-2way, silicone coated

2way standard:12Fr-24Fr,5-10ml / 30ml balloon

2way standard:26Fr-30Fr,30ml balloon

2way pediatric:6Fr-10Fr,3-5ml balloon

Available with the option of plastic valve or rubber valve


Packing: sterile individual blister

**Model:HD-DIS026**

Latex Foley-1 way, silicone coated

Size: 8Fr-24Fr

Packing: sterile individual blister

**Model: HD-DIS027**

All silicone foley catheter

2 way pediatric, 8Fr-10Fr, 3ml balloon

2 way standard, 12Fr-14Fr, 5-10/5-15ml balloon

2 way standard, 16Fr-24Fr, 5-10/15/30ml balloon

3 way irrigation(standard), 16Fr-24Fr, 30ml balloon

Packing: sterile individual blister

**Model: HD-DIS029**

Tracheotomy tube with cuff

Size: I.D.5.0mm-10.0mm

Packing: sterile individual blister

**Model: HD-DIS030**

Tracheotomy tube without cuff

Size: I.D.3.0mm–10.0mm

Packing: sterile individual blister

**Model: HD-DIS031**

Oral Pharyngeal Airway, Guedel style

Made from semirigid non-toxic Polyethylene

color-coded bite back

Guedel style enclosed center channel

Size: 40mm, 50mm, 60mm, 70mm, 80mm, 90mm,

100mm, 110mm, 120mm

Packing: sterile individual

**Model: HD-DIS032**

Oral Pharyngeal Airway, Berman type

Made from semirigid non-toxic colour coded for size identification

Size: 40mm,50mm,60mm,70mm,80mm,90mm,100mm,110mm

Packing: sterile individual

**Model: HD-DIS033**

Nasopharyngeal Airway

Size: I.D.12Fr-34Fr


Packing: sterile individual

**Model: HD-DIS034**

External male condom

Size: S, M, L, XL

Packing: sterile individual polybag OR blister

**Model: HD-DIS035**

Simple oxygen mask

With adjustable nose clip, elastic strap, 2.1m tubing

Size: XL,L,M,S

Packing: sterile individual polybag

**Model: HD-DIS036**

Venturi mask

In seven level: 24% 26% 28% 30% 35% 40% 50%

With adjustable nose clip, elastic strap, 2.1m tubing

Size: XL, L, M, S

Packing: sterile individual polybag

**Model: HD-DIS037**

Oxygen mask with reservoir bag

With adjustable nose clip, elastic strap, 2.1m tubing

Size: XL, L, M, S

Packing: sterile individual polybag

**Model: HD-DIS038**

Nebulizer with mouthpiece

2.1m corrugated tube

Packing: sterile individual polybag

**Model: HD-DIS039**

Nasal oxygen cannula

With 2.1m tube

Size: S, M, L

Packing: sterile individual polybag

**Model: HD-DIS040**

Nebulizer with mask

100% atomized granules (less than or equal to 5um)

With 2.1m tube

Size: XL, L, M, S

Packing: sterile individual polybag

**Model: HD-DIS041**

Connecting tube with yankauer handle

Option: Crow tip, Standard tip

Option: with control vent, without control vent

With connecting tube

Packing: sterile individual blister



Model: HD-DIS042

Scalpel

Size: 10, 11, 12, 12B, 13, 14, 15, 15C, 16,
18, 19, 20, 21, 22, 22A, 23, 24, 25, 36

Material: carbon steel / stainless steel

Packing: sterile individual




Model: HD-DIS043

Surgical blade

Size: 10, 10A, 11, 12, 12B, 13, 14, 15, 15C
16, 18, 19, 20, 21, 22, 22A, 23, 24, 25, 36

Material: carbon steel, stainless steel

Packing: sterile individual



HD-DIS044

Model: HD-DIS044

Skin graft blade

Multifunctional in surgical medical use

Sterile



HD-DIS045

Model: HD-DIS045

Double edge and single blade razor

Sterile


Packing: 50pcs/box

**Model: HD-DIS046**

Single edge and single blade razor

Sterile


Packing: 100pcs/box

**Model: HD-DIS047**

Single edge and double blade razor

Sterile

Packing: 100pcs/box



Needle Shape: straight, 1/4circle, 3/8circle or curved, J Needle and compound curved.

METRIC	mm	CATGUT (USP)	SILK (USP)	NYLON (USP)	POLYESTER (USP)
0.1	0.010–0.019			11/0	
0.2	0.020–0.029			10/0	
0.3	0.030–0.039			9/0	
0.4	0.040–0.049		8/0	8/0	
0.5	0.050–0.069		7/0	7/0	
0.7	0.070–0.099		6/0	6/0	6/0
1	0.100–0.149	6/0	5/0	5/0	5/0
1.5	0.150–0.199	5/0	4/0	4/0	4/0
2	0.200–0.244	4/0	3/0	3/0	3/0
3	0.300–0.349	3/0	2/0	2/0	2/0
3.5	0.350–0.399	2/0	0	0	0
4	0.400–0.499	0	1	1	1
5	0.500–0.599	1	2	2	2
6	0.600–0.699	2	3	3	
7	0.700–0.799	3	4	4	
8	0.800–0.899	4	5	5	
9	0.900–0.999	5	6		

Model: HD-DIS048

Suture



Model: HD-DIS049

Blood Lancet, Twist Top plastic handle

Model: HD-DIS050

Blood Lancet, Flake stainless steel



Model: HD-DIS051

Safety lancet

Needle: 18G, 21G, 26G

Different color available

Lancing Depth: 1.8mm, 2.4mm

Packing: 100pcs/box

**Model: HD-DIS052**

Lancet Device

Monitor diabetes patients glucose levels

Easily and safely use with an adjustable dial

Adjustable lancing depth

Any standard lancet can be used with these devices

Packing: individual case

**Model: HD-DIS053**

No adjustable safety lancet

Needle: 18G, 21G, 26G

Different color available

Lancing Depth: 1.8mm, 2.2mm, 2.8mm

Packing: 100pcs/box

**Model: HD-DIS054**

Adjustable safety lancet

26G needle, blue color

Lancing depth adjustable: 1.2mm, 1.8mm, 2.4mm

Packing: 100pcs/box

**Model: HD-DIS055**

Wooden tongue depressor

Adult use & child use

Size: 150mm X 18mm X 1.6mm

150mm X 19mm X 1.6mm

140mm X 14mm X 1.6mm

Non-sterilized

Packing: 100pcs/box

**Model: HD-DIS056**

single packed wooden tongue depressor

Single packed, sterile

Adult use

Size: 150mm X 18mm X 1.6mm

150mm X 19mm X 1.6mm,

140mm X 14mm X 1.6mm

Packing: sterile single-faced transparent OR double-faced paper

**Model: HD-DIS057**

Plastic Tongue depressor

Adult use

Size: 150mm X 19mm X 2mm

Non-sterilized

Packing: 100pcs/box

**Model: HD-DIS058**

Cervical scraper

With round tip & with tongue-shape tip

Size: 178mm x 18mm x 1.6mm

Non-sterilized or sterilized by ethylene oxide

Packing: single-faced transparent OR double-faced paper

**Model: HD-DIS059**

Applicator

Wooden-wand

Without cotton tip & with cotton tip

Big-head, little-head is available

Single-head with cotton tip & double-head with cotton tip

Size: 3", 6"

Non-sterilized

**Model: HD-DIS060**

Applicator sterile

Wooden-wand

Without cotton tip & with cotton tip

Big-head, little-head is available

Single-head with cotton tip & double-head with cotton tip

Size: 3", 6"

sterilized by ethylene oxide

Packing: sterile single-faced transparent OR double-faced paper

**Model: HD-DIS061**

Urine Bag

1500ml, 2000ml

Without bottom valve

**Model: HD-DIS062**

Urine Bag, 2000ml
With screw valve

**Model: HD-DIS063**

Urine Bag, 2000ml
With pull-push valve

**Model: HD-DIS064**

Urine Bag, 2000ml
With T valve

**Model: HD-DIS065**

Leg Bag
Capacity: 600ml, 750ml, 800ml, 1000ml



HD-DIS066

Model: HD-DIS066

Pediatric Urine Collector

Capacity: 100ml, 200ml



HD-DIS067

Model: HD-DIS067

Disposable vaginal speculum

Size: L, M, S


Packing: polybag / film / bulk

E.O. gas sterile

Type available: side screw, middle screw, French type,

South Asian type,

Spanish type, Australia type, American type




Model: HD-DIS075

Single-use spinal needle

Size: # 4 1/2-#12(26G-18G) Length 50-150mm

Use: Used in subarachnoid puncture for spinal anesthesia



Model: HD-DIS076

Dental Needle

Steriled by EO

Non-toxic, pyrogen free

Safety and convenience

It consists of adapter, cannula cap etc.

All sizes and all types are available

**Model: HD-DIS077**

Patient ID Bracelet

Adult/ infant

Color: pink, blue, white etc.

**Model: HD-DIS078**

Umbilical cord clamp

Made of non-toxic, medical-grade ABS/PE

Available colors: white, blue etc.

Packing: sterile individual polybag OR blister



Model: HD-DIS079

Pregnancy Test Strip (HCG)



Model: HD-DIS080


Pregnancy Test Cassette (HCG)



HD-DIS081

Model: HD-DIS081

Pregnancy Test Midstream (HCG)



HD-DIS082

Model: HD-DIS082

Drugs of Abuse Tests

–Detect AMP, BAR, BUP, BZO, COC, KET, MET, MDMA, MOP, MTD, OPI, OXY, PCP, PPX, TCA or THC

- High sensitivity
- Simple, easy and accurate
- Built-in control system
- Can read results at 5 minutes

1) Strip Format

–Lowest cost test format

2) Device Format

–Most commonly used test format

3) Panel Format

Convenient format for multi-drug screening

4) Multi-Drug Device Format

–May be any combination of 2–12 drug tests

5) Cup

–Most convenient format for urine collection and testing

–May be any combination of 2–12 drug tests

Rapid Test Name/Format/Specimen	Code	Sensitivity
AMP One Step Amphetamine Test Strip (Urine)	HD-DAM-101	1000 ng/ml
AMP One Step Amphetamine Test Device (Urine)	HD-DAM-102	1000 ng/ml
AMP300 One Step Amphetamine Test Strip (Urine)	HD-DAM-A101	300 ng/ml
AMP300 One Step Amphetamine Test Device (Urine)	HD-DAM-A102	300 ng/ml
BAR One Step Barbiturates Test Strip (Urine)	HD-DBA-101	300 ng/ml
BAR One Step Barbiturates Test Device (Urine)	HD-DBA-102	300 ng/ml
BZO One Step Benzodiazepines Test Strip (Urine)	HD-DBZ-101	300 ng/ml
BZO One Step Benzodiazepines Test Device (Urine)	HD-DBZ-102	300 ng/ml
COC One Step Cocaine Test Strip (Urine)	HD-DCO-101	300 ng/ml
COC One Step Cocaine Test Device (Urine)	HD-DCO-102	300 ng/ml
COC150 One Step Cocaine Test Strip (Urine)	HD-DCO-U101	150 ng/ml
COC150 One Step Cocaine Test Device (Urine)	HD-DCO-U102	150 ng/ml
KET One Step Ketamine Test Device (Urine)	HD-DKE-102	1000 ng/ml
MET One Step Methamphetamine Test Strip (Urine)	HD-DME-101	1000 ng/ml
MET One Step Methamphetamine Test Device (Urine)	HD-DME-102	1000 ng/ml
MET300 One Step Methamphetamine Test Strip (Urine)	HD-DME-A101	300 ng/ml
MET300 One Step Methamphetamine Test Device (Urine)	HD-DME-A102	300 ng/ml
MET500 One Step Methamphetamine Test Strip (Urine)	HD-DME-U101	500 ng/ml
MET500 One Step Methamphetamine Test Device (Urine)	HD-DME-U102	500 ng/ml
MDMA One Step Ecstasy Test Strip (Urine)	HD-DMD-101	500 ng/ml
MDMA One Step Ecstasy Test Device (Urine)	HD-DMD-102	500 ng/ml
MOP One Step Morphine Test Strip (Urine)	HD-DMO-101	300 ng/ml
MOP One Step Morphine Test Device (Urine)	HD-DMO-102	300 ng/ml
MTD One Step Methadone Test Strip (Urine)	HD-DMT-101	300 ng/ml
MTD One Step Methadone Test Device (Urine)	HD-DMT-102	300 ng/ml
OPI One Step Opiate Test Strip (Urine)	HD-DOP-101	2000 ng/ml
OPI One Step Opiate Test Device (Urine)	HD-DOP-102	2000 ng/ml
OXY One Step Oxycodone Test Strip (Urine)	HD-DOX-101	100 ng/ml
OXY One Step Oxycodone Test Device (Urine)	HD-DOX-102	100 ng/ml
PCP One Step Phencyclidine Test Strip (Urine)	HD-DPC-101	25 ng/ml
PCP One Step Phencyclidine Test Device (Urine)	HD-DPC-102	25 ng/ml
TCA One Step Tricyclic Antidepressants Test Strip (Urine)	HD-DTC-101	1000 ng/ml
TCA One Step Tricyclic Antidepressants Test Device (Urine)	HD-DTC-102	1000 ng/ml
THC One Step Marijuana Test Strip (Urine)	HD-DTH-101	50 ng/ml
THC One Step Marijuana Test Device (Urine)	HD-DTH-102	50 ng/ml
Multi-Drug Screen Test Panel (Urine)	HD-DOA-124	N/A
Multi-Drug Screen Test Cup (Urine)	HD-UCP-001	N/A



Model: HD-DIS083

Infectious Diseases Tests

- Aid in the detection of HBsAg, HBsAb, HBeAb, HBeAg, HBcAb, HBV, HCV, HIV, H. pylori, Malaria P.f., Strep A, Mononucleosis, TB, Dengue, Gonorrhoea, Chlamydia, Syphilis, Rotavirus, Adenovirus and ToRCH
- High sensitivity and specificity
- Simple, easy and accurate
- Built-in control system
- Can read results in 5-15 minutes

1) Strip Format

- Low cost test format

2) Device Format

- Most commonly used test format

3) Multi-test panel

- Best choice for multiple screening

Rapid Test Name/Format/Specimen	Code	Sensitivity
Chlamydia Rapid Test Device (Swab/Urine)	HD-ICH-502S	See Insert
HBV One Step Hepatitis B Virus Combo Test Device (Serum/Plasma)	HD-IHB-355	See Insert
HBsAb One Step Hepatitis B Surface Antibody Test Strip (Serum/Plasma)	HD-IHBsb-301	10 mIU/mL
HBsAb One Step Hepatitis B Surface Antibody Test Device (Serum/Plasma)	HD-IHBsb-302	10 mIU/mL
HBsAg One Step Hepatitis B Surface Antigen Test Strip (Serum/Plasma)	HD-IHBsg-301	1 ng/mL
HBsAg One Step Hepatitis B Surface Antigen Test Device (Serum/Plasma)	HD-IHBsg-302	1 ng/mL
HBsAg One Step Hepatitis B Surface Antigen Test Strip (Whole blood/Serum/Plasma)	HD-IHBsg-w301	0.5 ng/mL
HBeAb One Step Hepatitis B Envelope Antibody Test Device (Serum/Plasma)	HD-IHBeb-302	See Insert
HBeAg One Step Hepatitis B Envelope Antigen Test Device (Serum/Plasma)	HD-IHBeg-302	See Insert
HBcAb One Step Hepatitis B Core Antibody Test Device (Serum/Plasma)	HD-IHBcb-302	See Insert
HCV One Step Hepatitis C Virus Test Strip (Serum/Plasma)	HD-IHC-301	See Insert
HCV Hepatitis C Virus Rapid Test Device (Serum/Plasma)	HD-IHC-302	See Insert
HCV Hepatitis C Virus Rapid Test Strip (Whole Blood/Serum/Plasma)	HD-IHC-401	See Insert
HCV Hepatitis C Virus Rapid Test Device (Whole Blood/Serum/Plasma)	HD-IHC-402	See Insert
HIV Rapid Test Strip (Whole Blood/Serum/Plasma)	HD-IHI-401	See Insert
HIV Rapid Test Device (Whole Blood/Serum/Plasma)	HD-IHI-402	See Insert
H. pylori One Step Test Device (Serum/Plasma)	HD-IHP-302	See Insert
One Step H. pylori Antigen Test Device (Feces)	HD-IHP-602	See Insert
Malaria P.f. Rapid Test Strip (Whole Blood)	HD-IMA-401	See Insert
Malaria P.f. Rapid Test Device (Whole Blood)	HD-IMA-402	See Insert
Syphilis One Step Test Device (Serum/Plasma)	HD-ISY-U302	See Insert
Syphilis Rapid Test Strip (Serum/Plasma)	HD-ISY-U401	See Insert
Syphilis Rapid Test Strip (Whole Blood/Serum/Plasma)	HD-ISY-U411	See Insert
Syphilis Rapid Test Device (Whole Blood/Serum/Plasma)	HD-ISY-U402	See Insert
Dengue IgG/IgM Rapid Test Device (Whole Blood/Serum/Plasma)	HD-IDE-402	See Insert
Dengue IgG/IgM Rapid Test Device (Whole Blood/Serum/Plasma)	HD-IDE-401	See Insert
TB Tuberculosis Rapid Test Strip (Whole Blood/Serum/Plasma)	HD-ITB-401	See Insert
TB Tuberculosis Rapid Test Device (Whole Blood/Serum/Plasma)	HD-ITB-402	See Insert



HD-DIS084

Model: HD-DIS084**HIV-1/2 Test**

HIV-1/2 test is a immunochromatographic test for the qualitative detection of all isotypes antibodies specific to HIV-1 and HIV-2 simultaneously, in human serum, plasma or whole blood.

Features:

- The 3rd Generation Method (Direct Sandwich Method, Ag-Ab-Ag)
- Serum, Plasma, Whole Blood specimen
- Highly sensitive even during early infection stage
- Detect HIV-1 and HIV-2
- Sensitivity : 100%, Specificity : 99.8%
- Long shelf life : 24 months at room temperature

Composition:

- HIV-1/2 test device / strip/ POCT
- Assay diluent
- Lancet, alcohol swab and capillary pipette(20 μ /drop) optional



HD-DIS085

Model: HD-DIS085

Hepatitis B surface antigen (HBsAg) test

HBsAg test is designed as an initial qualitative test to detect Hepatitis B surface antigen (HBsAg) in human serum/plasma/whole blood at a concentration equal to or higher than 1 ng/ml. Hepatitis B is one of the major causes of hepatitis in humans, usually transmitted by close personal contact. Individuals carrying this virus may not be sick, but a source of infection. The HBsAg test is for healthcare professional use only and is available in cassette and dip strip formats.

HBsAg cassette and dip strip test	
Method	Lateral flow test
Sensitivity	> 1 ng/ml
Sample	Serum, plasma, whole blood
Read time	15–30 minutes
Storage	Room temperature
Shelf life	24 months from date of manufacture
Feature	Benefit
Easy, one-step procedure	Quick results; minimizes user error
Detects HBsAg at a concentration of higher than or equal to 1 ng/ml	Highly sensitive
Small sample volume required	Sample can be run reliably with minimal amount of specimen
24-month shelf life at room temp. (from date of manufacture)	No special storage requirements to maintain stability
Built-in Control	Indicates the test has been performed correctly
Made in China	Reliable, on-time supply & ready technical support



Model: HD-DIS089

Latex Obstetric Gloves

Anatomic shape for relaxed fit and prevention of fatigue
 Roughened surface on palm & finger tips for better grip
 Long & beaded cuff to reinforce cuff area, to maintain sterility
 Usage for Obstetric, Gynecology, Laboratory, Veterinary, etc

Sterilized by Gamma

Size: S/6.5", M/7.5", L/8.5"

Color: Natural

Length: 550mm

Thickness: 0.28mm (single surface)

Standard: ASTM D 3577, EN455, AQL=1.5

Packing: Sterile dual paper packing,

25 pairs/box, 4 boxes/case

**Model: HD-DIS090**

Disposable Penlight light for Diagnostic
Packing: individual box

**Model: HD-DIS091**

Disposable anesthesia breathing circuit with breathing
bag Hoses, length about 60''

Mask Elbow with Sampling Port

Fixed wye Piece

Bacterial Filter (Removable)

2.0/3.0L Breathing Bag, Natural Rubber Latex

Ventilator Hose

Different color of tubes and bags are available