



Model: HD-DRS001

Bandage Tubular Support

Material: 82% cotton+10% elastic rib+8% terylene

GSM: 530-580 g/m²

Tubigrip A, 4.50cmx10m

Tubigrip B, 6.25cmx10m

Tubigrip C, 6.75cmx10m

Tubigrip D, 7.50cmx10m

Tubigrip E, 8.75cmx10m and etc.

**Model: HD-DRS002**

Bandage Elastic Net Tubular

Material: cotton

Size: 0.5#, 1#, 2#, 3#, 4#, 5#, 5.5#, 5.8#, 6#, 7#, 8#

**Model: HD-DRS003**

Natural orthopedic padding

Material: fabric

Width: 5cm, 7.5cm, 10cm, 15cm, 20cm etc.

Length: 2.7m, 3m, 4.6m etc.

**Model: HD-DRS004**

Synthetic orthopedic padding

Material: sponge

Width: 5cm, 7.5cm, 10cm, 15cm, 20cm etc.

Length: 2.7m, 3m, 4.6m etc.

**Model: HD-DRS005**

Plaster of Paris bandage

Width: 5cm, 7.5cm, 10cm, 15cm, 20cm etc.

Length: 2.7m, 3m, 4.6m etc.

**Model: HD-DRS006/007**

Confirming Bandage

Color: bleach/unbleached

Width: 3cm, 4cm, 5cm, 7cm, 10cm etc.

Length: 4.5m etc.

**Model: HD-DRS008**

Elastic Bandage

Color: flesh, unbleached

Width: 2.5cm, 5cm, 7.5cm, 10cm etc.

Length: 4.5m etc.



HD-DRS010

Model: HD-DRS010

High elastic force bandage

Composition: 10% latex silk & 90% nylon

Color: flesh

Width: 2.5cm, 5cm, 7.5cm, 10cm etc.

Length: 4.5m etc.



HD-DRS011

Model: HD-DRS011

Cotton Crepe bandage

Color: flesh, bleach or unbleached

Width: 2.5cm, 5cm, 7.5cm, 10cm etc.

Length: 4.5m etc.

**Model: HD-DRS012**

Non-woven self-adhesive bandage

Width: 2.5cm, 5cm, 7.5cm, 10cm etc.

Length: 4.5m etc.

**Model: HD-DRS013**


All cotton self-adhesive bandage

Width: 2.5cm, 5cm, 7.5cm, 10cm etc.

Length: 4.5m etc.

**Model: HD-DRS013**

Woolen Bandage

**Model: HD-DRS015**

Zinc oxide plaster

Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm etc.

Length: 1m, 5m, 10m etc.

Color: white, flesh

Packing: simplified packing OR Plastic cylinder OR Tinplate shell

**Model: HD-DRS016**


Microporous surgical tape non-woven

Width: 1.25cm, 2.5cm, 5cm, 7.5cm etc.

Length: 5Y, 10Y etc.

Available Packing:

- Simplified packing
- Plastic cylinder
- Spool Blister
- With cutting device (Dispenser)

**Model: HD-DRS017**


Transparent surgical tape PE

Width: 1.25cm, 2.5cm, 5cm, 7.5cm etc.

Length: 5Y, 10Y etc.

Available Packing:

- simplified packing
- Plastic cylinder
- Spool Blister
- With cutting device (Dispenser)



Model: HD-DRS018

Wound adhesive plaster

Option:

Flexible Cotton wound adhesive plaster

Colored and carton wound adhesive plaster

Non-woven wound adhesive plaster

PU PE PVC wound adhesive plaster

Wound adhesive plaster with special shape



Model: HD-DRS019

Strip Cap

Size: 18", 21", 24" etc.

Color: green, blue, white and red etc.

GSM: 12g etc.

**Model: HD-DRS020**

Doctor Cap

Type option: tie type, elastic type

Size: 18", 21", 24" etc.

Color: green, blue, white etc

GSM: 20G, 30G etc.

**Model: HD-DRS021**

Bouffant Cap

Size: 18", 21", 24" etc.

Color: green, blue, white etc.

GSM: 12G etc.

**Model: HD-DRS022**

N95 Face Mask

Option: with valve, without valve

**Model: HD-DRS023**

Dust-proof face mask

Prevents spread of germs and filters out dust and pollen

Size: 160g, 180g, 200g

**Model: HD-DRS024**

Paper Mask

Type: single ply, two ply

GSM: 14,16G etc.

Size: 20*7cm

**Model: HD-DRS025**

Non-woven surgical face mask

Type: earloop or tie

2ply or 3ply is available

Color: green, blue, white and pink etc

GSM: 43G (2ply), 63G (3ply) etc.

**Model: HD-DRS026**

Face Mask with Eye Shield

Options: 3 ply 63g (pp18-mb20-pp25)+PC shield, earloop type

4 ply 103g (pp18-activated carbon

40-mb20-pp25), tie type

Shield shape option: rectangle-shaped, shield size 18x13cm

Fan-shaped, shield size 30x12.5cm

**Model: HD-DRS027**

Shoe cover

Material: CPE/PE, 2.5g/pc, 2.7g/pc etc. size 15x36cm etc.

Material: Non-woven, GSM 30g, 35g, 40g etc,
size 14 x40cm etc.

**Model: HD-DRS028**

Non-skid Shoe Cover

Size: 36cmx15cm, 38cmx16cm, 40cmx17cm etc

GSM:35G, etc.

**Model:HD-DRS029**

Surgeon gown

Material: SMS, spunlacer

Size: M, L, XL

Color: blue, white etc.

Sterile & non-sterile

GSM: 20G, 25G, 40G etc.

**Model: HD-DRS030**

Isolation gown

Material: SMS, PP

Size: M, L, XL

Color: blue, white etc.

GSM: 30G, 40G, 50G etc.

**Model: HD-DRS031**

Over coat

Material: SMS, PP

Size: M, L, XL

Color: blue, white etc.

GSM: 30-50G etc.



HD-DRS032

Model: HD-DRS032

Visit Gown

Material: PP

Color: blue, white, green and yellow

GSM: 40G etc



HD-DRS033

Model: HD-DRS033

Gauze Swab

40's x 40's cotton yarn

With or without X-Ray thread

Mesh: 19x15, 24x20 etc.

Size: 5cmx5cm, 7.5cmx7.5cm, 10cmx10cm, 10cmx20cm etc.

Ply: 8, 12, 16 etc.

Color: white, green

Sterile or non sterile

Fold or unfold

Paper-paper packing OR paper-plastic packing

**Model: HD-DRS034**

Gauze Roll

40's x 40's cotton yarn

With or without X-Ray thread

Mesh: 19 x 15, 24 x 20 etc.

Size: 36" x 50M, 4Ply; 36" x 100Y, 4Ply etc.

**Model: HD-DRS035**

Wow gauze bandage

40's x 40's cotton

Mesh: 19 x 15, 24 x 20 etc.

Size: 5cm x 4cm, 7.5cm x 4cm, 10cm x 4cm etc.

Individual cellophane packing

**Model: HD-DRS036**

Triangle bandage

100% cotton or rayon fabric cotton

Yarn: 40's x 40's

Mesh: 40x34, 50x30, 48x48 etc.

Size: 36" x 36" x 51", 40" x 40" x 56"

With or without softy pin

Colored: grey or bleached

**Model: HD-DRS037**

Gauze swab (single packed)

40's x 40's cotton yarn

With or without X-Ray thread, sterile

Mesh: 19x15, 24x20 etc.

Size: 5cm x 5cm, 7.5cm x 7.5cm, 10cm x 10cm,

10cm x 20cm etc.

Ply: 8, 12, 16 etc.

Color: white, green



Model: HD-DRS038

Alcohol prep pad

70% Isopropyl alcohol non-woven

Size: 65mm x 30mm, 65 x 56mm, 90x90mm etc.



Model: HD-DRS039

Laparotomy sponge

40's x 40's cotton yarn

With or without X-Ray thread

With or without cotton loop

Mesh: 19 x 15, 24 x 20 etc.

Size: 30cm x 30cm, 45cm x 45cm etc.

Ply: 4, 8, 12 etc

**Model: HD-DRS040**

Absorbent cotton wool

Weight: 50g, 100g, 200g, 500g etc.

**Model: HD-DRS041**

PE APRON

Thickness: 0.018mm, 0.020mm etc.

Option: embossed or smooth

Size: 700 x 1150mm, 810 x 1500mm, etc.

Color: white, blue etc.

**HD-DRS042**

Vinyl Glove (P.V.C Glove)

Color: Clear, blue, green, red etc.

Powdered & Powder free

Size: XS, S, M, L, XL

**Model: HD-DRS043**

Sterile Latex Surgical Gloves

Powered or powered-free

Size: 5.5#, 6#, 6.5#, 7#, 7.5#, 8#, 8.5#, 9#

**Model: HD-DRS044**

PE Glove

Material: LDPE, HDPE etc.

Option: embossed or smooth

Size: S, M, L, XL

Thickness: 0.01–0.035mm

Color: white, blue etc.

**Model: HD-DRS045**

PE Glove with long sleeve

Length: 40–93cm

Thickness: 0.02–0.04mm

Option: embossed or smooth