



Model:HD-HYP001

Disposable syringe luer slip

1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 50ml

3 parts, non-toxic, pyrogen free

Packing: sterile individual polybag OR blister

**Model: HD-HYP002**

Disposable syringe luer lock

2ml, 3ml, 5ml, 10ml, 20ml etc.

3 parts, non-toxic, pyrogen free

packing: sterile individual polybag OR blister

**Model: HD-HYP003**

Disposable insulin syringe

1ml, separated needle

Non-toxic, pyrogen free

packing: sterile individual polybag OR blister

**Model: HD-HYP004**

Disposable insulin syringe

1ml, with two jackets

with incoipoiated needle

Non-toxic, pyrogen free

Packing: sterile individual polybag OR blister

**Model: HD-HYP005**


Disposable insulin syringe

1ml, with two jackets

with separated needle

Non-toxic, pyrogen free

Packing: sterile individual polybag OR blister

**Model: HD-HYP006**

Disposable tuberculin syringe

1ml, with incoipoated needle

Non-toxic, pyrogen free

Packing: sterile individual polybag OR blister

**Model: HD-HYP007**

Disposable vaccination syringe

1ml, IU-100

Non-toxic, pyrogen free

Packing: sterile individual polybag OR blister

**Model: HD-HYP008**

Disposable syringe catheter tip

50/60ml, without needle

Non-toxic, pyrogen free

Packing: sterile individual polybag OR blister

**Model: HD-HYP009**

Disposable hypodermic needle

Needle size : 19-27G, different length accordingly

Non-toxic, pyrogen free

Packing: sterile individual polybag OR blister

**Model: HD-HYP010**

Safety auto-destroy syringe, disposable
1ml , 3ml , 5ml , 10ml , 30ml

Non-toxic, pyrogen free

Luer slip OR luer lock (except 1ml)

Packing: sterile individual polybag OR blister

**Model: HD-HYP011**

Disposable infusion set

Vented spike with cap, small chamber with filter


PE or ABS adjustor

PVC tube length 150cm

Rubber luer lock or luer slip

With needle, Non-toxic, pyrogen free

Packing: sterile individual polybag OR complex film

**Model: HD-HYP012**

Disposable infusion set

Vented spike with cap, large chamber with filter

PE or ABS adjustor

PVC tube length 150cm

With needle, Non-toxic, pyrogen free

Packing: sterile individual polybag OR complex film

**Model: HD-HYP013**

Disposable infusion set

Vented spike with cap, large chamber with filter

Y-injection set

PE or ABS adjustor

PVC tube length 180cm

Rubber luer lock with cap

With needle, Non-toxic, pyrogen free

Packing: sterile individual polybag OR complex film

**Model: HD-HYP014**

Disposable blood transfusion set

Spike with cap, without outlet

Chamber with filter

PE or ABS adjustor

PVC tube 150cm

With needle, Non-toxic, pyrogen free

Packing: sterile individual polybag OR complex film

**Model: HD-HYP015**

Disposable blood transfusion set

Vented spike with cap, chamber with filter

PE or ABS adjustor

PVC tube 150cm

Spike with cap, with outlet

With needle, Non-toxic, pyrogen free

Packing: sterile individual polybag OR complex film

**Model: HD-HYP016**

Disposable burette set

For children use

100ml, 150ml

Vented spike with valve

PE or ABS adjustor, cylinder chamber

PVC tube, Y-injection set

Non-toxic, pyrogen free

Packing: sterile individual polybag

**Model: HD-HYP017**

Disposable scalp vein set

International colour code


Butterfly needle, 18-27G luer slip

25cm length tube

Matched with infusion set or transfusion set

Non-toxic, pyrogen free

Packing: sterile individual polybag OR blister



Model: HD-HYP018

Disposable scalp vein set
 International colour code
 Butterfly needle 18-27G, luer lock
 25cm length tube
 Matched with infusion set or transfusion set
 Non-toxic, pyrogen free
 Packing: sterile individual polybag OR blister



Model: HD-HYP019

I.V. Cannula
 Option: Pen-like model
 With wings model
 With injection valve model

Color	Gauge	O.D.x Length:(mm)
Yellow	24G	0.7x19
blue	22G	0.8x25
Pink	20G	1.0x33
Green	18G	1.2x45
Grey	16G	1.6x55
Red	14G	2.0x55

Packing: sterile individual Polybag OR blister packing



Model: HD-HYP020

Heparin Cap/Stopper
 Closing stopper for intermittent injections
 Silicon rubber
 At least 20 times injection guaranteed
 Luer lock fitting
 Packing: sterile individual blister packing

Model: HD-HYP021

Heparin cap
 Closing stopper for intermittent injection
 Silicon rubber
 At least 20 times injection guaranteed
 Transparent/visible
 Luer lock fitting
 Packing: sterile individual blister packing



Model: HD-HYP022

3-way stopcock
 Two female luer port and one male taper port
 The cap clicks as it stops precisely at each position
 to guarantee accuracy turning 360°
 The cap can be colour coded in blue or red
 Non-toxic, pyrogen free
 Packing: sterile individual blister



HD-HYP023

Model: HD-HYP023

2-way stopcock


Two female luer port and one male taper port

The cap clicks as it stops precisely at each position to guarantee accuracy turning 180°

The tap can be colour coded in blue or red

Non-toxic, pyrogen free

Packing: sterile individual blister



HD-HYP024

Model: HD-HYP024

3-way stopcock with extension set

Pre-assembled extension tube, OD: 2.6mm

Length: 100cm, 200cm, 300cm, 500cm, etc

For infusion management or for pressure monitor

Non-toxic, pyrogen free

Packing: sterile individual blister

**Model: HD-HYP025**

Extension Tube

Made of clear, Non-toxic PVC, OD: 2.6mm

Length: 100cm, 200cm, 300cm, 500cm, etc

For infusion management, for pressure monitor

Non-toxic, pyrogen free

Packing: sterile individual blister