



Model: HD-XRP001

Lead Rubber Jacket/Lead Vest/Lead Apron Set
/Lead Coat Apron

Size: 1000 × 600mm, 900 × 600mm

Front 0.35 mmPb or 0.50mmPb options

Sleeve and Back 0.25mmPb



Model: HD-XRP002

Hands Protective



Model: HD-XRP003

Arms Protective



Model: HD-XRP004

Lead Gloves



HD-XRP005

Model: HD-XRP005

Intervient Radiation Protective Gloves
(0.025-0.03mmPb)



HD-XRP006

Model: HD-XRP006

Protective Collar



HD-XRP007

Model: HD-XRP007

Protective Collar



HD-XRP008

Model: HD-XRP008

Protective Collar



HD-XRP009

Model: HD-XRP009


Protective Collar



HD-XRP010

Model: HD-XRP010

Profiled Collar




Model: HD-XRP011

Universal X-ray Protective Glasses



Model: HD-XRP012

Red Eyepiece



Model: HD-XRP013

Protect Scrotum



Model: HD-XRP014
Protective Short Skirt



Model: HD-XRP015
Protective Triangle



Model: HD-XRP016
Protective Cap




Model: HD-XRP017

Combined Protective Scarf Apron



Model: HD-XRP018

Lead Rubber



Model: HD-XRP019

Cassette Cover Plate



Model: HD-XRP020

Mobile Swivel Arm Apron Rack



Model: HD-XRP021

Chest Thick Ruler



Model: HD-XRP022

Film Viewer



HD-XRP023

Model: HD-XRP023

Film Viewer



HD-XRP024

Model: HD-XRP024

X-Ray Film, Green/Blue



HD-XRP025

Model: HD-XRP025

Lead Screen



Model: HD-XRP026

Lead Chair



Model: HD-XRP027

Lead Films-storing Box



Model: HD-XRP028

Wall/Standing type X-ray Film Cassette Shelf



Model: HD-XRP029

Standing X-ray Film Cassette Shelf



Model: HD-XRP030

X-Ray Envelope